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Epithesen-
versiegelung matt

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Instructions for use

Matt sealing agent for
epitheses

Before using this product, please read these instructions for use carefully!

The manufacturer will not accept any liability for damage resulting from non-compliance with these instructions for use.

1. Description of the product	
Matt sealing agent for epitheses is a surface sealing agent for silicone surfaces. Additionally, integrated UV-protection avoids fading of the silicone surface.	
2. Indication	
Stains on surfaces, such as Multisil stains, are fixed by matt sealing agent for epitheses and protected against wearing off.	
3. Contraindication	
Based on the current state of knowledge there are no known contraindications. The product should not be used for patients with a suspected allergy until allergy tests have been completed to confirm that the patient is not allergic to the materials contained in the product.	
4. Safety and hazard information	

Matt sealing agent for epitheses – primer	
Signal word:	Danger
Hazard statements:	
H225	Highly flammable liquid and vapor.
H260	In contact with water releases flammable gases which may ignite spontaneously.
H304	May be fatal if swallowed and enters airways.
H332	Harmful if inhaled.
Precautionary statements:	
Prevention:	
P210	Keep away from heat, sparks, hot surfaces and other ignition sources. No smoking.
P233	Keep container tightly closed.
P261	Avoid breathing vapor.
Response:	
P301+P310	IF SWALLOWED: Immediately call a POISON CENTER or doctor / physician.
P303+P361+P353	IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower.
P331	Do NOT induce vomiting.
P335+P334	Brush off loose particles from skin. Immerse in cool water.

Product contains: Solvent naphtha (petroleum), light aliphatic	
Matt sealing agent for epitheses – Coating	
Signal word:	Danger
Hazard statements:	
H225	Highly flammable liquid and vapor.
H319	Causes serious eye irritation.
Precautionary statements:	
Prevention:	
P210	Keep away from heat, sparks, hot surfaces and other ignition sources. No smoking.
P233	Keep container tightly closed.

Response:	
P303+P361+P353	IF ON SKIN (or hair): Remove / Take off immediately all contaminated clothing. Rinse skin with water / shower.
P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P337+P313	If eye irritation persists: Get medical advice / attention.

5. Storage and durability	
Matt sealing agent for epitheses – Primer and Coating:	
Storage:	
P403+P235	Keep cool and store in a well-ventilated place.

6. Processing	
6.1 Use	
The primer of matt sealing agent for epitheses is applied thinly onto the epithesis surface using a brush. Then the epithesis is characterized with Multisil stains. A uniform layer of coating of matt sealing agent for epitheses is dabbed on carefully. Multisil Epithetik fibers can be applied onto the coating layer, which is still soft. The coating material of matt sealing agent for epitheses is dried (e.g. with a hairdryer) for two minutes at 65 °C. The matting powder is applied after 15 minutes (avoid skin contact) and placed into hot water (70 to 80 °C) for two minutes. Matt sealing agent for epitheses increases the adhesion within 24 hours.	

7. Symbols	
7.1 Symbols according to DIN EN ISO 15223-1	
REF	Reference number
LOT	Lot number
Can be used until YYYY-MM	
Upper temperature limit: 28 °C	
Keep away from sunlight	
Manufacturer	

7.2 Symbols according to GHS regulation	
Flamme	GHS02, flame
!	GHS07, exclamation mark
GHS08, health hazard	

7.3 Symbols according to NFPA (USA)	
1 0 HEALTH FIRE REACTIVITY	Hazard diamond

8. Additional information	
The product may only be used by dentists, dental technicians and adequately trained dental staff. Only original tools and components may be used for processing. The lot numbers of all components should be recorded to ensure traceability and to handle complaints. These instructions for use are based on state-of-the-art methods and equipment and our own experience. The product may only be used for the indication described under item 2. The user himself is responsible for processing the product. Liability for incorrect results shall be excluded since the manufacturer does not have any influence on further processing. Any occurring claims for damages may only be made up to the value of our products.	

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